

SUMMARY

Reduction in technical repeat and recall rate after implementation of artificial intelligence driven quality improvement software

Eby PR, Martis L, Paluch J, Jang J, Chan A. 2022 RSNA Annual Meeting, Chicago, USA.

Why it matters

This large study affirms that unnecessary imaging can be reduced with use of mammography quality improvement software. Volpara Analytics™ is effective in improving quality and reducing repeats and recalls for large breast cancer screening centers that produce a high-volume of mammograms.

Key takeaways

This study is the largest North America–based mammographic image quality (IQ) evaluation to date. Following installation of Volpara Analytics™ software, significant improvements of 6% and 8% were observed in objectively measured breast positioning and compression IQ, respectively. This corresponded to a concomitant 78% reduction in technical repeats and recalls (TR), from 0.77% to 0.17%.



Study location

Virginia Mason Franciscan Health, Seattle, WA



Study dataset

>48K exams; 42 technologists across 9 clinics



Study design

Quality improvement evaluation of IQ indicators and TR rates following installation of Analytics software

Methods

- Analytics software was installed in 2019 to objectively measure breast positioning and compression IQ for all images and all technologists.
- Technologists reviewed interactive dashboards that provided continuous IQ feedback, trend analyses, and benchmarking, enabling them to identify focus areas for improvement and set goals
- TR rates (n=40 technologists; 48,143 exams) and IQ indicators (n=42 technologists; 48,874 exams) were compared for the first (“Baseline”, April 2019 - March 2020) versus the most recent (“Current”, April 2021 - March 2022) 12-month periods following Analytics installation.

Study results

- Comparing *Baseline* versus *Current* periods:
 - TR rates reduced by 78% (0.77% to 0.17%)
 - Mean overall Quality Score increased by 6% (2.28 to 2.42)
 - Proportion of images meeting Target Compression increased by 8% (59.1% to 63.6%)
 - Proportion of images scored Perfect or Good increased by 6% (56.4% to 59.8%)
- Analytics software facilitated mammographic IQ evaluation on an unprecedented scale and has the potential to improve outcomes for both providers and consumers of mammography screening.