

# Imaging System Compatibility

## VolparaDensity Releases 3.3.1, 3.3.2, 3.3.3, 3.4, 3.4.1, and 3.4.2

Volpara Health's products are compatible with a broad range of mammography systems for the main providers of full-field digital mammography, digital breast tomosynthesis, and computed radiography systems.

For more information on DICOM requirements for each imaging system, see the Volpara Solutions DICOM Conformance Statement, available for download from our Support page.

Manufacturer	Model	3D	2D	CR
Canon	MGU-1000D		>	
Fujifilm	AMULET Innovality/ASPIRE Cristalle	>	>	
	AMULET/ASPIRE HD		>	
	AMULET s/ASPIRE HDs		>	
	AMULET f/ASPIRE HD Plus		>	
	ClearView CSm			>
	FCR PROPECT CS Plus			>
	Innovality/Cristalle	>	>	
GE	Senographe Pristina	>	>	
	SenoClaire	>		
	Senographe Essential   DS   2000D		>	
Hologic	Selenia Dimensions	>	>	
	Selenia		>	

Manufacturer	Model	3D	2D	CR
IMS Giotto	Giotto Image SDL   3DL		>	
Medi-Future	Brestige		>	
	SOUL		>	
Metaltronica	Helianthus		>	
Phillips	MicroDose L50   L30		>	
	MammoDiagnost DR		>	
	Mammo Diagnost		>	
PlanMed	Nuance		>	
	Nuance Excel		>	
Sectra	MDM 1.5		>	
Siemens	Mammomat Fusion		>	
	Mammomat Inspiration	>	>	
	Mammomat Inspiration PRIME		>	
	Mammomat NovationDR		>	
	Mammomat Revelation	>	>	

# Volpara Scorecard

## Release 4.1

Manufacturer	Model	3D	2D	CR
Canon	MGU-1000D		>	
Fujifilm	AMULET/ASPIRE HD		>	
	AMULET s/ASPIRE HDs		>	
	AMULET f/ASPIRE HD Plus		>	
	ClearView CSm			>
	FCR PROPECT CS Plus			>
GE	Senographe Pristina	>	>	
	SenoClaire	>		
	Senographe Essential   DS   2000D		>	
Hologic	Selenia Dimensions	>	>	
	Selenia		>	
IMS Giotto	Giotto Image SDL   3DL		>	
Medi-Future	SOUL		>	
Metaltronica	Helianthus		>	
Phillips	MicroDose L50		>	

Manufacturer	Model	3D	2D	CR
PlanMed	Nuance		>	
	Nuance Excel		>	
Sectra	MDM 1.5		>	
Siemens	Mammomat Fusion		>	
	Mammomat Inspiration	>	>	
	Mammomat Inspiration PRIME		>	
	Mammomat NovationDR		>	
	Mammomat Revelation	>	>	

# End of Life (EOL) Announcements

Manufacturer	Model	Product Type	Final Mfg Date	Final Mfg Date	"End of Life" Plan	Available Replacement Product
Matakina	VolparaAnalytics (all versions)	DataAnalytics	Jan 2015	Jan 2015	<p>Technical Support available until Feb 2024.</p> <p>Service Replacement parts are limited to available inventory. No available upgrade path.</p>	VolparaAnalytics
Volpara Health Technologies	VolparaAnalytics (all versions)	DataAnalytics	Mar 2017	Oct 2017	<p>Technical Support available until Feb 2024.</p> <p>Service Replacement parts are limited to available inventory. No available upgrade path.</p>	VolparaAnalytics